



MAK and Isolation

The MAK equipment will follow the JDH Infection Control process. Located in section 1.11 cleaning, Sterilization and Disinfection of Patient Care Equipment

MAK Hints

- ❖ **Position the cart within your line of sight-but outside of the room if possible**
- ❖ **If you bring the cart into isolation room wipe down with JDH approve cleaning materials**
- ❖ **Place scanner in clear plastic bag prior to entering the isolation**
- ❖ **Remove scanner from plastic bag before leaving room**

Do not override a patient on isolation